





		Par	Parts List		
Item Qty	y Description 2	Description	Material	Vendor	Note
	INNER SS TUBE	3/4 OD STAINLESS STEEL TUBE - 0.083 WALL	- 0.083 WALL STAINLESS STEEL		
2	OUTER AL PIPE	AL 6061 SCH 40 PIPE	AI 6061 PIPE		ANODIZE BLUE
3 1	BOTTOM CUP		ALUMINUM 6061		ANODIZE BLUE
4	TEFLON RING	MACHINED MECHANICAL GRADE PTFE	TEFLON		
<u>ე</u>	TEFLON HANDLE	MACHINED MECHANICAL GRADE PTFE	TEFLON		
6 2	MODIFED DOWEL	MODIFIED $1/8 \times 3/8$ DOWEL PIN	STAINLESS STEEL	MSC 67600122 or EQUIV.	MSC 67600122 or EQUIV. PRESSED & WELDED TO SS TUBE
7 1	DOWEL PIN	PIN — GROUND DOWEL, 1/4 x 1 1/4	STAINLESS STEEL	MSC 67601047 or EQUIV.	BRIGHT FINISH, CHAMFER OK
8 2	O-RING	AS 568 - 018	TEFLON	MSC 31955313 or EQUIV.	EQUIV. TIGHTER TOLERANCES THAN AS568A
9	SPRING	COMPRESSION .975 x .072 x 1.5	STAINLESS STEEL	MSC 03309218 or EQUIV.	
10 1	ENGRAVED LABEL	0.056" THK. ROWMARK RED/WHITE	ACRYLIC		LASERMARK ENGRAVING STOCK
11 2	SCREW $\#4-40 \times 3/8$	HEXAGON SOCKET HEAD CAP SCREW	STAINLESS STEEL	MSC 05664040 or EQUIV.	
12 2	SCREW $\#4-40 \times 5/8$	HEXAGON SOCKET HEAD CAP SCREW	STAINLESS STEEL	MSC 05664065 or EQUIV.	
13 3	SCREW $\#2-56 \times 1/4$	SLOTTED PAN HEAD MACHINE SCREW	STAINLESS STEEL	MSC 67568725 or EQUIV.	


				ISSUED	APPR.	CHECK	DRAWN 07.30.2003 Mitch Miller	CONTRACT NO -	FILE NAME Transfer-Handel-RevB.dwg
SCALE 1:1	t	U I	SIZE FSCM NO		PARTS LIST		)	ICSOI/	でいる
WEIGHT			DWG NO				.	ICSOI / SIMIB / SSIKIL	/ では、この / でであった。
SHEET 3 OF 14	t		REV			TANULE	]	L/SLAIC	Ì












Notes:


- 1. RAW MATERIAL 0.056" THK. ROWMARK ACRYLIC ENGRAVING STOCK. RED BACKGROUND WITH WHITE LETTERS.
- 2. TEXT SAYS "PUSH DOWN TO TURN" WITH AN ARROW INDICATING THE DIRECTION FOR "OFF" AND "ON". A 3-DIGIT SERIAL NUMBER IS ENGRAVED IN A SMALLER FONT UNDER "TO TURN".
- 3. LARGER TEXT HEIGHT IS TO BE 11±0.5 POINTS (0.153±0.007")ENGRAVED IN A SAN SERIF BOLD FONT (E.G. ARIAL BLACK OR ARIAL NARROW BOLD).
- 4. THE SERIAL NUMBER TEXT HEIGHT IS TO BE  $7.5\pm0.5$  POINTS (0.104 $\pm0.007$ ") ENGRAVED USING A SAN SERIF FONT OF THE SAME TYPE FAMILY (E.G. ARIAL OR ARIAL NARROW).
- 5. LAYOUT APPROVAL IS REQUIRED IF USING DIFFERENT FONTS FROM THOSE LISTED ABOVE.
- 6. AN ADOBE ILLUSTRATOR OR AUTOCAD .DWG FILE CAN BE SUPPLIED WITH THE TEXT IF REQUESTED (PLEASE SPECIFY VERSION).
- 7. ALWAYS CONFIRM THE BEGINNING AND ENDING SERIAL NUMBERS BEFORE STARTING THE RUN.
- 8. ASSEMBLY: SEE SHEET 12




,XX ±0.010	ALL DIMENSIONS ARE IN INCHES	IF IN DOUBT, ASK, TITTED	THIRD ANGLE	<b>\</b>				
			ISSUED	APPR.	CHECK	DRAWN 07.30.2003 Mitch Miller	CONTRACT NO -	FILE NAME Transfer-Handel-RevBawg
SCALE 2:1	AM_B -	SIZE FSCM NO				)	JICSO /	
WEIGHT	I	DWG NO		D LABEL			CSU/SIMIB/SSIKIL/SILAIC	
SHEET 9 🛮 F 14	B	REV			HANDLE	- - - - - - 1	/ SLAIC	)

- Notes: 1. RAW MATERIAL MECHANICAL GRADE PTFE TEFLON
- 2. I.D. TOLERANCE IS TIGHTER THAN AS-568A SPEC.
- Ÿ. SPECS FOR AS-568A-018 O-RINGS I.D. NOMINAL 3/4" I.D. 0.739±0.009 C.S. 0.070±0.003
- C.S. O.D. NOMINAL NOMINL 1/16" 7/8"
- PART# 31955313). O-RINGS IN A PACKAGE OF AS-568A-018 TEFLON O-RINGS (MSC 4. WE FOUND THAT WE ONLY REJECTED APPROX. 40-50 % OF THE
- 5. WE TESTED THE O-RINGS ON A GO/NO-GO TESTER. O-RINGS WERE ACCEPTED IF THE MOVED FREELY OVER A  $\emptyset$ 0.741 OD TUBE.
- 6. IF THE ID OF THE 0-RINGS IS TOO SMALL, THE HANDEL WILL SEIZE WHEN USED UNDER LIQUID NITROGEN WHEN THE TEFLON CONTRACTS AGAINST THE WALLS OF THE STAINLESS STEEL TUBE.
- 7. WE MADE OUR TESTER FROM A SHORT PIECE OF 3/4" STAINLESS STEEL TUBING WHICH HAD BEEN TURNED DOWN. MOST OF THE TUBE WAS TURNED DOWN TO \$\pi\_0.735\$. THE OTHER END WAS LEFT FULL SIZE. THIS ALLOWED US TO QUICKLY TEST O-RINGS. IF THEY MOVED EASILY ACROSS THE \$\pi\_0.741\$ SECTION, THEY COULD BE RETAINED ON THE TUBE BY THE FULL DIAMETER SECTION. ONES THAT DID NOT MOVE FREELY OVER THE \$\pi\_0.741\$ SECTION WERE PLACED IN A SEPARATE BIN.
- ASSEMBLY: SEE SHEET 14

œ




±0.010	ALL DIMENSIONS ARE IN INCHES	IF IN DOUBT ASK.	THIRD ANGLE	<b>&gt;</b>				
			ISSUED	APPR.	CHECK	DRAWN 07.30.2003 Mitch Miller	CONTRACT NO -	FILE NAME Transfer-Handel-RevB.dwg
SCALE 5:1	_AM_B —	SIZE FSCM NO			CAVVE		ILSU,	
WEIGHT		DWG ND		□-RING -018			CSUI/SIMIBI/SSIKIL/SILAIC	
SHEET				$\infty$				7 / 7
SHEET 10 OF 14	В	REV			HANJLE			



ASSEMBLY NOTES: ASSEMBLY OF THE ENGRAVED LABEL INTO THE TEFLON HANDLE IS ASFOLLOWS:

- 1. THE ENGRAVED LABEL (10) IS TO BE EMBEDDED INTO THE TOP OF THE TEFLON HANDLE (5) WITH LOCTITE HYSOL E-30CL ADHESIVE (50 ML CARTRIDGE IS LOCTITE #29329). THIS EPOXY BASED ADHESIVE HAS A 30 MINUTE POT LIFE.
- 2. AIR MUST BE PURGED FROM MIXING NOZZLE BEFORE BEGINNING ASSEMBLY. CARE MUST ME USED WHEN APPLYING THE EPOXY TO MINIMIZE THE INTRODUCTION OF AIR BUBBLES.
- 3. FIRST APPLY A SMALL AMOUNT OF ADHESIVE INTO THE TOP OF THE TEFLON HANDLE (5). THIS IS TO FILL THE BLIND HOLES AND UNDER THE LIP OF THE DOVETAIL.
- 4. NEXT PRESS THE ENGRAVED LABEL INTO RECESS IN THE TEFLON HANDLE. WORK AROUND SLIGHTLY TO ALLOW EPOXY TO COME UP OVER THE EDGES OF THE LABEL AND FILL DOVETAIL.
- 5. THE INSERT IS TO BE LEFT CLOSE TO CENTERED IN THE RECESSED OPENING WITH THE ENGRAVED LABEL ROTATED SUCH THAT FLAT IN THE HANDLE IS BELOW THE TEXT AND THE #40-40 CLEARANCE HOLES IN THE TEFLON HANDLE ARE PERPENDICULAR TO THE TEXT DIRECTION FROM ABOVE AND PARALLEL TO THE TEXT DIRECTION TO THE RIGHT AND THE BLIND HOLES TO THE BOTTOM AND LEFT (AS SHOWN).
- 6. NOW FILL EPOXY OVER THE TOP OF THE LABEL UNTIL A SLIGHT POSITIVE MENISCUS IS FORMED. CARE SHOULD BE TAKEN TO AVOID THE INTRODUCTION OF AIR BUBBLES. A COUPLE (0-3) SMALL AIR BUBBLES MAY BE UNAVOIDABLE, BUT LARGE NUMBERS OR LARGE SIZE AIR BUBBLES ARE NOT ACCEPTABLE.
- 7. SET ASIDE TO CURE. NOTE THAT THE POT LIFE IS 30 MINUTES BUT IT TAKES OVER 24 HOURS FOR IT TO CURE. ALSO WE FOUND THAT THE SURFACE TAKES LONGER TO CURE FULLY. THE CLARITY OF THE SURFACE CAN BE DAMAGED BY FINGERPRINTS, OILS, OR PRESSURE FOR 3—7 DAYS. THE PARTS CAN BE HANDLED BEFORE THIS TIME WITH CARE NOT TO TOUCH THE EPOXY SURFACE. WE TYPICALLY MAKE AN EPOXY TEST PIECE TO DETERMINE WHEN THE PARTS HAVE CURED SUFFICIENTLY FOR HANDELING. CURING CAN BE ENHANCED IN LOW HUMIDITY ENVIRONMENTS AND AT ELEVATED TEMPERATURES (NOT TO EXCEED 176°F / 80°C).
- 8. MAKE SURE EPOXY IS CURED BEFORE COMPLETING ASSEMBLEY STEPS SHOWN ON SHEET 14.
- 9. IF YOU HAVE ANY QUESTIONS ABOUT THE PROCESS, PLEASE ASK. WE CAN LOAN YOU A SAMPLE TO DEMONSTRATE THE FINISHED EPOXY EMBEDDING IF NECESSARY.






## JICSG / SIMIB / SSIRIL / SILAIC

CASSETTE TRANSFER HANDLE

EXPLODED VIEW - HANDLE LABEL

SCALE 2:1	AM_B -	SIZE FSCM NO
WEIGHT	1	DWG NO
SHEET		
12 OF		
14	₩	REV



ASSEMBLY NOTES - OUTER TUBE:

- 1. PRESS TEFLON RING (4) ONTO BOTTOM CUP (3) AND ATTACH WITH  $\#2-56 \times 1/4$ " SLOTTED PAN HEAD MACHINE SCREW (13).
- 2. SLIDE BOTTOM CUP (3) ONTO OUTER ALUMINUM PIPE (2) AND ATTACH WITH  $\#4-40 \times 3/8$ " SOCKET HEAD CAP SCREWS (11).

ILE NAME Transfer-Handel-RevBdwg  INCSG / SIMB / SSIRIL / SILAIC  INTRACT NO -  INTRAC	SHEET 13 OF	WEIGHT	SCALE 1:2	
RACT NO -  107.30.2003 Mitch Miller  CA:  SIZE FS  R -			ţ	
RACT NO - CA			<del>B</del> 1	
RACT NO -  1 07.30.2003 Mitch Miller  CA		DWG NO	SIZE FSCM NO	
RACT NO - CA:				ISSUED
NAME Transfer-Handel-RevB.dwg  RACT NO -  07.30.2003 Mitch Miller  CA:	ASSEME	VIEW - DUTER	EXPLODED	APPR.
	HANI		LAVVEII	CHECK
9	- >		) ) )     	DRAWN 07.30.2003 Mitch Miller
9	/ SLA	SMIB / SSIKIL	IICSO /	CONTRACT NO -
	/ ~ " "			FILE NAME Transfer-Handel-RevB.dwg

 $\bowtie$ 

